

# Accuracy counts for all oral appliances.



**ProGauge**  
AUTOCLAVEABLE

#### **FIT AND MONITOR ANY ORAL APPLIANCE.**

The ProGauge device assists in accurately taking a bite registration for all mandibular advancement appliances. Use it to record your patient's protrusive range of motion for titration purposes both clinically or in a sleep lab. The ProGauge from Airway Management ensures proper hardware placement. It also allows a precise vertical setting from 2mm to 12mm.



866.264.7667  
TAPINTOSLEEP.COM