


# TAP® 3 Elite

## Thornton Adjustable Positioner® with double bar



**The TAP® 3 Elite double bar is indicated specifically for your Medicare patients.**

Features include:

- a fixed mechanical hinge and inseparable pivot (double bar/hook)
- in-built mechanism that allows for adjustment in increments of ¼ mm
- can protrude the mandible beyond the front teeth to maximum protrusion and farther to maximum passive protrusion
- easily adjusted by the beneficiary while in the mouth, using a key
- retains the adjustment setting when removed
- maintains treatment position during sleep and prevents dislodgement of the device

*“Your Patients hold the key”*



**For more information or to place an order, call Customer Service at 866.264.7667**

