

TAP®

Thornton Adjustable Positioner® *with closed hook*



The TAP® closed hook is indicated specifically for your Medicare patients.

Features include:

- a fixed mechanical hinge and inseparable pivot (closed hook/lingual bar)
- in-built mechanism that allows for adjustment in increments of ¼ mm
- can protrude the mandible beyond the front teeth to maximum protrusion and farther to maximum passive protrusion
- easily adjusted by the beneficiary with a convenient single point adjustment knob, while in the mouth
- retains its adjustment setting when removed
- maintains treatment position during sleep and prevents dislodgement of the device

For more information or to place an order, call Customer Service at 866.264.7667

